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### Conditions for Figures 1-3

column: Ascentis Express Silica, 10 cm x 2.1 mm I.D.,  
2.7 µm particles (53930-U)  
mobile phase: 10:90, (A) 100 mM ammonium formate, pH 3.0 with concentrated formic acid (B) acetonitrile  
temperature: 35 °C  
flow rate: 0.4 mL/min.  
detection: UV at 254 nm  
injection volume: 1 µL  
1. Acenaphthene, 80 mg/mL in mobile phase  
2. Adenosine, 35 mg/mL in mobile phase  
3. Cytosine, 75 mg/mL in mobile phase

Figure 1. Typical QA Chromatogram for Ascentis Express HILIC Column


Figure 2. Chromatogram from First Injection on a New Column (wrong technique)



Figure 3. Chromatogram After Switching Wash Solvent Bottles Between Weak and Strong Solvents

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## Reference

1. Waters On-Line Help: Selecting an Acquity Injection Technique.

© 2012 Sigma-Aldrich Co. LLC. All rights reserved. SIGMA, SAFC, SIGMA-ALDRICH, ALDRICH, and SUPELCO are trademarks of Sigma-Aldrich Co. LLC, registered in the US and other countries. FLUKA is a trademark of Sigma-Aldrich GmbH, registered in the US and other countries. SAFC brand products are sold through Sigma-Aldrich, Inc. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see product information on the Sigma-Aldrich website at [sigma-aldrich.com](http://sigma-aldrich.com). Ascentis is a registered trademark of Sigma-Aldrich and Sigma-Aldrich Biotechnology, LP. Fused-Core is a registered trademark of Advanced Materials Technology, Inc. Acquity UPLC and Waters are registered trademarks of Waters.

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